



BIOREADER® 6000-F-C

Model [Part no.]	Purpose	Fluoro filter-changer	Camera Mpixel	DELL PC 3 years NBD on-site service and warranty	DELL Monitor 3 years replacement warranty
Pro-F-C b [2507]	Elispot/Zones/ Plaque/CFU Colony counting. Typicall in 96-24 well plates and petridishes (35mm- 100 mmm).	Capability of a color colony counter 4 filter blocks (i.e. 3 fluoro filters +1 enzymatic filter). Brightfield, toplight and darkfield. UV excitation for fluorogenic colonies and other assays.	10	I7 processor, 4 GB RAM, 1 TB GB hard drive RAID 0, NAS server	24" 1920* 1600 Pixel
Pro-F-C g [2507]	Same like Pro-F-C b	6 filter blocks (5 filters maximal +1 enzymatic automatic operated)	10	I7 processor, 4 GB RAM, 1 TB GB hard drive RAID 0, NAS server	24" 1920* 1600 Pixel