



BIOREADER® 6000-F

Model [Part no.]	Purpose	Fluoro-filter changer	Camera Mpixel	DELL PC 3 years NBD on-site service and warranty	DELL Monitor 3 years replacement warranty
Pro-F a [1735]	Routine	3 filter blocks (i.e. 2 Fluoro filters +1 enzymatic filter) * Automatic operated changer	1.2	Dual core, 4 GB RAM, 320 GB hard drive	24'' 1920 * 1200 Pixel
Pro-F b [1584]	Research + routine	4 filter blocks (i.e. 3 Fluoro filters +1 enzymatic filter) * Automatic operated changer	3	Dual core, 4 GB RAM, 500 GB hard drive	24'' 1920 * 1200 Pixel
Pro-F g [1366]	21 CFR part 11 Clinical study, research	6 filter blocks (5 filters maximal +1 enzymatic Automatic revolver operated)	5	Quad core, 4 GB RAM, 1 TB hard drive, RAID 0, NAS	24'' 1920 * 1200 Pixel