



BIOREADER® 6000-V

Model [Part no.]	Purpose	User interface	Camera Mpixel	DELL PC 3 years NBD on-site service and warranty	DELL Monitor 3 years replacement warranty
Pro-V a [2221]	To be used for Plaque/CFU/Elispot: in 96 and 12/24 or 96 and 6 well plates , or 96 well and 48 well plates, or 12/24 and 6 well plates , or 48 and 6 well plates, or 48 and 12/24 well plates (should be specified)	'Ease of use' 'professional'	1.2	Dual core, 4 GB RAM, 160 GB hard drive	20" 1600 * 1200 Pixel
Pro-V b [1674]	Elispot/Plaque/CFU/MLA/Micro dilution: in 96 ,24/12 and 6 well plates , or 96 , 48, 12/24 well plates (should be specified)	'Ease of use', 'professional' 'GLP'	3	Dual core, 4 GB RAM, 320 GB hard drive	20" 1600 * 1200 Pixel
Pro-V g [1733]	Elispot/Plaque/CFU/MLA/Micro dilution/Colonies/ Antibiotic potency in 10 cm petridishes: <ul style="list-style-type: none"> • 96 and 12/24 and 6 and petridishes, • or 96 and 48 and 12/24 and 6 well plates 	'Ease of use', 'professional', 'GLP'	5	Quad core, 4 GB RAM, 500 GB hard drive, RAID 0, NAS	24" 1920 * 1200 Pixel
Optional	Larger case, if 100 mm petridish are to be used				