



BIOREADER® 6000-Z

Model [Part no.]	Purpose	User interface	Zoom/lens changer	Camera Mpixel	DELL PC 3 years NBD on-site service and warranty	DELL Monitor 3 years replacement warranty
Pro-Z a [1471]	Elispot, Plaque + CFU: in 96, 48, 24, 12, 6 well plates and 35- 100 mm petridishes. New: with Top light.	'Ease of use'	Software controlled zoom and focus	1	Dual core, 4 GB RAM, 160 GB hard drive	20" 1600 * 1200 Pixel
Pro-Z b [2222]	Elispot, Plaque + CFU: in 96, 48, 24, 12, 6 well plates and any petridish. New: with Top light.	'Ease of use', 'professional'	High definition motorized zoom. Top/bottom homogenous illumination	3	Dual core, 4 GB RAM, 320 GB hard drive	24" 1920 * 1200 Pixel
Pro-Z g [2211]	Elispot, Plaque, CFU, Infections biology, Genetic Toxicology, Antibiotic potency: in 96 and 12/24 (or 6/12) well plates. Optional: 35 mm - 10 cm petridishes [Item 2219].	'Ease of use', 'professional', 'GLP'	High definition motorized zoom. Top/bottom homogenous illumination	5	Quad core, 4 GB RAM, 500 GB hard drive, RAID 0, NAS	24" 1920 * 1200 Pixel